



ITEMS AVAILABLE FOR LOAN

(up to 3 months)

MOBILITY

KNEE SCOOTER

(Knee walker)



Used by individuals that can't bear weight on one of their lower ankles or feet.

WALKER WITH SEAT



Walker with brakes, seat, steerable wheels and can be folded up.

UPRIGHT WALKER



Supports the weight on the forearms instead of wrists and hands. This reduces pressure and pain in the lower back, neck, arms and wrists.

BARIATRIC WHEELCHAIR



For larger consumers or individuals with external fixators on hips or legs. Larger weight capacity, up to 500 pounds. 24-26” seat width, verify doors are wide enough for access.

STANDARD WHEELCHAIR



Standard seat width, 18-20 inch. Typically 250 pound weight capacity.

TRANSFER/TRANSPORT WHEELCHAIR



For attendant to transport consumer, not able to be self-propelled.

STAIR CANE



To assist when going up or down stairs.

4 PRONG CANE



Cane that has increased stability and can handle greater weight.

CANE



June 2024

CHAIR CANE



A chair cane offers a sturdy cane for walking and easily folds open for a quick place to sit and rest.

ICE PICK ATTACHMENT FOR CANES



Attaches to the end of a single point cane for traction.

CRUTCHES



3 WHEELED TRIANGLE WALKER



3 wheeled walkers are lighter weight, more maneuverable, and will fit into corners (like kitchen cupboards) closer as compared to a 4 wheeled walker. They do not have a seat and are less stable.

HEMI WALKER



One sided walker for those with hemiparesis.

STANDARD AND JUNIOR WALKERS



Both have several height adjustment settings.

STAND ASSIST HANDLE



Increase leverage to assist with standing, reduce risk of falls during transfers.

LEG LIFTER



Used to move or reposition leg after hip or knee surgery while in bed, car, wheelchair, etc.

HOME SAFETY & ACCESSIBILITY

BED RAIL



Slides under the mattress to provide additional assistance to get in/out of bed.

LEVER EXTENDER



**Attaches to recliner handle to allow for neutral grip,
good for those with arthritis.**

TOILET SAFETY FRAME



Attaches under the mounting bolts of the toilet seat to add arms for leverage to sit and stand.

LARGE KEY TELEPHONE AMPLIFIED WITH PHOTO BUTTONS



For those with hearing impairment, cognitive impairment, or low vision. Amplified handset and ringer volume. Can insert pics onto white keys to call, bright red++++ light when ringing.

REMOTE CONTROL OUTLETS




**Can turn on 3 different outlets with one remote.
Reduces the need to bend over or turn on lamps with
hard-to-turn twist switches for those with limited
hand dexterity.**

WRITING GUIDE



Device slides easily along writing surface with only slight downward pressure needed. Helpful for those with limited hand coordination or weak grasp.

VIBRATING CHIME PAGER



Will detect when a door or window is open and chime and vibrate the receiver up to 100 ft away.

KEY TURNER




Assist people with arthritis by providing a handle for their keys.

AUTO OPENING SCISSORS



For individuals with limited hand strength and mobility.

4 FT WHEELCHAIR RAMP



Aluminum portable foldable ramp. Carries like a suitcase.

HEALTH/MEDICAL

TALKING ALARM CLOCK FOR MEDS



Will cue a person with poor memory or executive function skills to take their meds up to 4 times per day, alarm will continue until stopped.

KITCHEN AIDS

KNOB AID



Molds itself over oven and other knobs to help people with arthritis rotate.

J POPPER



Pulls the tab on a can (food or beverage) to open.

JAR OPENER



Helpful for those with arthritis or anyone with a weakened grip.

GRIP IT-GROCERY BAG HOLDERS



Holds several bags with one grip for the user.

BOTTLE CAP OPENER



Provides an oversize sticky grip to turn smaller caps on bottles.

PERSONAL CARE
BATHING
DRESSING

BATHTUB RAIL



Attaches to edge of bathtub to assist in getting in and out of tub.

BATH STOOL



Typical weight capacity 300 pounds.

BATH STOOL WITH BACK



Typical weight capacity 300 pounds.

BATH TRANSFER BENCH



Weight limit 300-400 pounds.

SUCTION GRAB BAR



Attaches to smooth wall surfaces only. Will not work with textured or tiled surfaces.

SOCK AID



Helpful for those who find it hard to reach their feet.

BUTTON HOOK



Helpful for those with arthritis, use of only one hand, weak grip and pinch strength. Pulls button through the hole with one hand.

DOFF AND DONNER



**To put socks or compression stockings on.
Helpful How-To videos available on YouTube.**

DRESSING AID



Pull up pants, pull sleeves through shirts, push off socks. Helpful for people with limited mobility.

LONG HANDLED SHOE HORN



HEARING LOSS

POCKET TALKER



Detachable Mic



Belt Clip



Tone Control




Volume Control



Sound amplifier for individuals with hearing loss; the user wears headphones. Produces a very clear sound.

VISION LOSS

MAGNIFYER SHEET



Magnifies a whole page.

DOME MAGNIFYER



To magnify lines of text for those with low vision.

LARGE PRINT HIGH CONTRAST KEYBOARD




Helpful for those with low vision and/or decreased visual attention to tasks.

HANDHELD MAGNIFYER



Non-lighted, handheld magnifier.

LARGE KEYS CALCULATOR



Helpful for those with low vision or hand tremors.

READING GUIDE STRIP



Assists people with reading difficulties to focus on one line of print at a time through the orange window.

VEHICLE MOBILITY

AUTO ASSIST HANDLE



Attaches to car door to add a handle for entry and exit.

CAR GLARE SHIELD



Clips onto car visor to extend the shaded area.

GAS CAP REMOVER



To aid in removing gas caps, arthritis, decreased hand strength.

HANDY BAR



Fits into a car door latch to add an additional handle to grab, also has a seat belt cutter.

RECREATION

CARD HOLDERS



GARDENING TOOL SET



Large handled garden tools.